



Mouse Anti-Isofenphos-methyl monoclonal antibody, clone IFM-2 (CABT-L7792)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Tested Cross-Reactivity: Isocarbophos 2ppb
Immunogen	Isofenphos-methyl with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	IFM-2
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	C14H22NO4PS; Isofenphos-methyl is a broad spectrum, organophosphorus insecticide, used to control soil pests in grain and oil crops such as click beetles, mole cricket and grub.
Keywords	Isofenphos-methyl; IFM