



# Mouse Anti-Myclobutanil monoclonal antibody, clone MCB (CABT-L7793)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Myclobutanil with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	MCB
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

## BACKGROUND

**Introduction**

C<sub>15</sub>H<sub>17</sub>ClN<sub>4</sub>; CAS Number: 88671-89-0.

Myclobutanil is a light yellow solid used as a fungicide. It is a triazole chemical used as a fungicide. It is a steroid demethylation inhibitor, specifically inhibiting ergosterol biosynthesis.

---

**Keywords**

Myclobutanil; MCB

---