



Mouse Anti-Triapenthenol monoclonal antibody, clone TPT (CABT-L7797)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Triapenthenol with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	TPT
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Triapenthenol is a herbicide. It is highly soluble in water. Little is known regarding its environmental fate but it is considered to be persistent in soils. It is not highly toxic to mammals. It shows a high level of toxicity to birds. There is little information on its toxicity to biodiversity.

Keywords

Triapenthenol; TPT
