



# Mouse Anti-Gibberellic acid monoclonal antibody, clone GA (CABT-L7799)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Gibberellic acid with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	GA
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

## BACKGROUND

<b>Introduction</b>	Gibberellic acids (Gibberellins) are naturally occurring plant hormones that are used as plant growth regulators to stimulate both cell division and elongation that affects leaves and stems.
<b>Keywords</b>	Gibberellic acid; GA