



Mouse Anti-Forchlorfenuron monoclonal antibody, clone FCF (CABT-L7800)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Forchlorfenuron with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	FCF
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction Forchlorfenuron is a plant growth regulator.

Keywords Forchlorfenuron; FCF
