



Mouse Anti-Dimethomorph monoclonal antibody, clone DMM (CABT-L7803)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Dimethomorph with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DMM
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Dimethomorph is a fungicide with systemic function. It is used for treating mildew and root rot caused by organisms such as Pythium and Phytophthora species.
Keywords	Dimethomorph; DMM