



# Mouse Anti-Dimethachlon monoclonal antibody, clone DMC (CABT-L7804)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Dimethachlon with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DMC
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

## BACKGROUND

**Introduction** Dimethachlon. CAS# 24096-53-5.

---

**Keywords** Dimethachlon; DMC

---