



Mouse Anti-Procymidone monoclonal antibody, clone PCM (CABT-L7805)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Procymidone with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PCM
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Procymidone is a pesticide. It is often used for killing unwanted ferns and nettles, and as a dicarboximide fungicide for killing fungi, for example as seed dressing, pre-harvest spray or post-harvest dip of lupins, grapes, stone fruit, strawberries.

Keywords

Procymidone; PCM
