



# Mouse Anti-Chlordimeform monoclonal antibody, clone CDF (CABT-L7808)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Chlordimeform with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	CDF
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

## BACKGROUND

**Introduction** Chlordimeform is an acaricide (pesticide) active mainly against motile forms of mites and ticks and against eggs and early instars of some Lepidoptera insects.

---

**Keywords** Chlordimeform; CDF

---