



Mouse Anti-Phorate monoclonal antibody, clone PR (CABT-L7809)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Tested Cross-Reactivity: Disulfoton 40ppb
Immunogen	Phorate with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	PR
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Phorate is an organophosphate used as an insecticide and acaricide.
Keywords	Phorate; PR