



# Mouse Anti-Diuron monoclonal antibody, clone DR (CABT-L7811)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	DCMU with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	DR
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

## BACKGROUND

**Introduction**

C<sub>9</sub>H<sub>10</sub>Cl<sub>2</sub>N<sub>2</sub>O.

Diuron is a white crystalline solid. It is a wettable powder. The primary hazard is the threat to the environment. Immediate steps should be taken to limit its spread to the environment. It can cause illness by inhalation, skin absorption and/or ingestion. It is used as a herbicide.

---

**Keywords**

Diuron; DR

---