



Mouse Anti-Abietic acid monoclonal antibody, clone ATA (CABT-L7819)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Abietic acid with carrier protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Clone	ATA
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

MF: C₂₀H₃₀O₂; CAS: 514-10-3

Abietic acid is an abietane diterpenoid that is abieta-7,13-diene substituted by a carboxy group at position 18. It has a role as a plant metabolite. It is an abietane diterpenoid and a monocarboxylic acid. It is a conjugate acid of an abietate.

Keywords

Abietic acid; ATA
