



Mouse Anti-GM2 Monoclonal Antibody, Clone NL2-27 (CABT-L7827)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	GM2 (NeuAc)
Target	GM2
Immunogen	GM2(NeuAc)+Monophosphoryl lipid A+ Treharose dimycolate emulsion
Isotype	IgM, κ
Source/Host	Mouse
Clone	NL2-27
Purification	The antibody is purified from mouse ascites fluid.
Conjugate	Unconjugated
Applications	TLC, FC, ELISA, IHC
Format	Liquid
Concentration	Lot specific
Buffer	PBS with 0.05% sodium azide
Preservative	0.05% sodium azide
Storage	Store at -20°C.
Ship	Wet ice