



Mouse Anti-GD1a Monoclonal Antibody, Clone HNS28 (CABT-L7828)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Ganglioside GD1a (IV3 NeuAc α , II3 NeuAc α -Gg4 Cer)
Target	GD1a
Immunogen	Purified ganglioside GD1a
Isotype	IgM, κ
Source/Host	Mouse
Clone	HNS28
Purification	The antibody is purified from mouse ascites fluid.
Conjugate	Unconjugated
Applications	TLC, FC, ELISA, IHC
Format	Liquid
Concentration	Lot specific
Buffer	PBS with 0.05% sodium azide
Preservative	0.05% sodium azide
Storage	Store at -20°C.
Ship	Wet ice