



Mouse Anti-GM3 Monoclonal Antibody, Clone HNS7 (CABT-L7831)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Ganglioside GM3
Target	GM3
Immunogen	Purified ganglioside GM3 adsorbed to Salmonella Minnesota
Isotype	IgM, κ
Source/Host	Mouse
Clone	HNS7
Purification	The antibody is purified from mouse ascites fluid.
Conjugate	Unconjugated
Applications	TLC, FC, ELISA, IHC
Format	Liquid
Concentration	Lot specific
Buffer	PBS with 0.05% sodium azide
Preservative	0.05% sodium azide
Storage	Store at -20 or -80°C. It is not recommended storing at 4°C. Avoid freeze/thaw cycles.
Ship	Dry ice