



Mouse Anti-Lewis X Monoclonal Antibody, Clone X84-41 (CABT-L7835)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Lewis X
Target	Lewis X
Immunogen	Salmonella Minnesota coated with Lewis X-active glycolipids
Isotype	IgM, κ
Source/Host	Mouse
Clone	X84-41
Purification	The antibody is purified from mouse ascites fluid.
Conjugate	Unconjugated
Applications	TLC, FC, ELISA, IHC, WB, IP
Format	Liquid
Concentration	Lot specific
Buffer	PBS with 0.05% sodium azide
Preservative	0.05% sodium azide
Storage	Store at -20°C.
Ship	Wet ice