



Mouse Anti-Lewis Y Monoclonal Antibody, Clone I29B (CABT-L7836)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Lewis Y antigen
Target	Lewis Y
Immunogen	Immunized with human gastric cancer cells (MKN74), and screened using cancer glycolipids
Isotype	IgG3, κ
Source/Host	Mouse
Clone	I29B
Purification	The antibody is purified from mouse ascites fluid.
Conjugate	Unconjugated
Applications	TLC, FC, ELISA, IHC, WB
Format	Liquid
Concentration	Lot specific
Buffer	PBS with 0.05% sodium azide
Preservative	0.05% sodium azide
Storage	Store at -20°C.
Ship	Wet ice