



# Mouse anti-MTPAP monoclonal antibody, clone 2E4 (CABT-L7846)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	The immunogen used to generate this antibody corresponds to human mtPAP.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	2E4
<b>Purification</b>	Affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC/IF, IHC-P Recommended Dilution WB 1:500-1:3000 ICC/IF 1:100-1:1000 IHC-P 1:100-1:1000
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	1XPBS, 20% Glycerol (pH7)
<b>Storage</b>	Store at -20°C.
<b>Ship</b>	Wet ice