



Mouse anti Shiga toxin type II monoclonal antibody, clone Mab55 (CABT-L7852)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Cap: Mouse anti Shiga toxin type II monoclonal antibody, clone Mab55
Specificity	Shiga Toxins II
Target	E. coli Shiga Toxins II
Immunogen	Purified Native Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Shiga toxin-producing E.coli
Clone	Mab55
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA(Cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L7852 - CABT-L7853
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide

Storage	Store at -20°C.
Ship	Wet ice