



Mouse anti Staphylococcus aureus SEA monoclonal antibody, clone Mab302 (CABT-L7854)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Cap: Mouse anti Staphylococcus aureus SEA monoclonal antibody, clone Mab302
Specificity	SEA specific
Immunogen	Purified Native Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	S. aureus
Clone	Mab302
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA(Cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L7854 - CABT-L7855
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide

Storage	Store at -20°C.
Ship	Wet ice