



# Mouse anti Staphylococcus aureus SEA monoclonal antibody, clone Mab275 (CABT-L7855)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Det: Mouse anti Staphylococcus aureus SEA monoclonal antibody, clone Mab275
<b>Specificity</b>	SEA specific
<b>Immunogen</b>	Purified Native Protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	S. aureus
<b>Clone</b>	Mab275
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(Det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L7854 - CABT-L7855
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% sodium azide

<b>Storage</b>	Store at -20°C.
<b>Ship</b>	Wet ice