



Mouse anti *Staphylococcus aureus* SEC monoclonal antibody, clone Mab314 (CABT-L7861)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Det: Mouse anti <i>Staphylococcus aureus</i> SEC monoclonal antibody, clone Mab314
Specificity	SED specific
Immunogen	Purified Native Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. aureus</i>
Clone	Mab314
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA(Det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L7860 - CABT-L7861
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide

Storage	Store at -20°C.
Ship	Wet ice