



Mouse anti Staphylococcus aureus SEE monoclonal antibody, clone Mab313 (CABT-L7862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Cap: Mouse anti Staphylococcus aureus SEE monoclonal antibody, clone Mab313
Specificity	SEE specific
Immunogen	Purified Native Protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	S. aureus
Clone	Mab313
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA(Cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L7862 - CABT-L7863
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide

Storage	Store at -20°C.
Ship	Wet ice