



Mouse Anti Norfloxacin monoclonal antibody, clone NFLX (CABT-L7864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Norfloxacin with carrier protein
Isotype	IgG
Source/Host	Mouse
Clone	NFLX
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in ELISA and LFIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction	Norfloxacin is used to treat certain types of infections, including infections of the urinary tract and prostate
Keywords	Norfloxacin; Noroxin; Chibroxin; Norfloxacin