



Mouse Anti Cefapirin monoclonal antibody, clone CFP (CABT-L7865)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CEP with carrier protein
Isotype	IgG
Source/Host	Mouse
Clone	CFP
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in ELISA and LFIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction

Cefalotin, cefapirin, cefazolin, cefadroxil and cefalexin are all first generation cephalosporins and are active against most aerobic Gram-positive cocci, many anaerobic cocci and bacilli, and some aerobic Gram-negative bacteria.

Keywords

Cefapirin; Cephapirin; Cephapirine; Cefaprin
