



# Mouse Anti Roxarsone monoclonal antibody, clone ROX (CABT-L7866)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Roxarsone with carrier protein
Isotype	IgG
Source/Host	Mouse
Clone	ROX
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in ELISA and LFIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store at -20°C.
Ship	Wet ice

## BACKGROUND

<b>Introduction</b>	Roxarsone is an organoarsenic compound that has been used in poultry production as a feed additive to increase weight gain and improve feed efficiency, and as a coccidiostat.
<b>Keywords</b>	Roxarsone; ROX