



Mouse Anti Virginiamycin monoclonal antibody, clone VGM (CABT-L7867)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Virginiamycin with carrier protein
Isotype	IgG
Source/Host	Mouse
Clone	VGM
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in ELISA and LFIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store at -20°C.
Ship	Wet ice

BACKGROUND

Introduction Virginiamycin is a streptogramin antibiotic similar to pristinamycin and quinupristin/dalfopristin. It is a combination of pristinamycin IIA (virginiamycin M1) and virginiamycin S1.

Keywords Virginiamycin; eskalinv; mikamycin; ostreogrycin; patricin; streptogramin
