



# Mouse Anti Nifurstyrenate monoclonal antibody, clone NFS (CABT-L7869)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Nifurstyrenate with carrier protein
<b>Isotype</b>	IgG
<b>Source/Host</b>	Mouse
<b>Clone</b>	NFS
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Suitable for use in ELISA and LFIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Purified, Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% sodium azide
<b>Storage</b>	Store at -20°C.
<b>Ship</b>	Wet ice

## BACKGROUND

**Keywords**Nifurstyrenate

---