



Rabbit Anti-Human Vimentin (p-Ser55) polyclonal antibody (CABT-L7873)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KLH conjugated synthesised phosphopeptide derived from human Vimentin around the phosphorylation site of Ser55: YA(p-S)SP
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in WB, FC, and ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at -20°C.
Ship	Wet ice