



Rabbit Anti-E. coli O149:K91,K88ac Polyclonal Antibody (CABT-L837R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Escherichia coli O149:K91,K88ac (Abbotstown)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E.coli
Purification	Unpurified
Conjugate	Unconjugated
Applications	Agg
Format	Liquid
Size	1 ml
Buffer	Antisera contains 0.01% Thiomersal.
Preservative	None
Storage	Store at 4°C.
Ship	Wet ice

BACKGROUND

Introduction	Escherichia coli (E. coli) are Gram negative enteric bacteria that live in the human gut. Its presence in nature is an indication of human fecal pollution. E. coli also causes human
---------------------	---

disease such as urinary tract infections and neonatal meningitis. Infection with particular enterotoxigenic strains can lead to life threatening intestinal diseases.

Keywords

Proteobacteria;Gammaproteobacteria;Enterobacteriales;Enterobacteriaceae;Escherichia;E. coli
