



Rabbit Anti-E. coli K99 polyclonal Antibody (CABT-L855R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	E.coli K99
Target	E. coli K99
Immunogen	Escherichia coli K99
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E.coli
Purification	Unpurified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	5 ml
Buffer	Antisera contains 0.01% Thiomersal.
Preservative	None
Storage	Store at 4°C.
Ship	Wet ice