



Mouse Anti-Human IL2RA monoclonal antibody, clone C-H4 (CABT-L875M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognises the Interleukin 2 Receptor a chain (IL-2 Ra)
Target	IL2RA
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	C-H4
Purification	Affinity purification
Conjugate	Unconjugated
Applications	ELISA Capture (Matched Pair) Recommended Working Concentration: ELISA: 0.5-2.5 µg/ml We recommend the following for sandwich ELISA (Capture - Detection): CABT-L875M - CABT-L876M; CABT-L875M - CABT-L876MU
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Preservative	None

Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.
----------------	--

Ship	Wet ice
-------------	---------

BACKGROUND

Keywords	IL2RA; interleukin 2 receptor; alpha; p55; CD25; IL2R; TCGFR; IDDM10; interleukin-2 receptor subunit alpha; TAC antigen; IL-2R subunit alpha; IL-2 receptor subunit alpha; HuMax-TAC
-----------------	--

GENE INFORMATION

Synonyms	IL2RA; interleukin 2 receptor; alpha; p55; CD25; IL2R; TCGFR; IDDM10; interleukin-2 receptor subunit alpha; TAC antigen; IL-2R subunit alpha; IL-2 receptor subunit alpha; HuMax-TAC
-----------------	--
