



Rabbit Anti-Mouse LY6A monoclonal antibody, clone KN64-22 (CABT-L933)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	LY6A
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Clone	KN64-22
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, IHC
Cellular Localization	Cell membrane.
Positive Control	Mouse brain tissue, mouse kidney tissue.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.
Preservative	0.05% Sodium Azide
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

BACKGROUND

Keywords

Ly-6A.2/Ly-6E.;Ly-6A.2/Ly-6E.1;Ly6a;LY6A_MOUSE;Ly6al;Lymphocyte antigen 6A-2/6E-1;Lymphocyte antigen Ly 6A.2/Ly 6E.1;SCA-1;Stem cell antigen 1;T cell activating protein;T-cell-activating protein;TAP antibody
