



Rabbit Anti-FGFR2 monoclonal antibody, clone KN21-71 (CABT-L942)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	FGFR2
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	KN21-71
Purification	Protein A purified.
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Cell membrane, Golgi apparatus, Cytoplasmic vesicle.
Positive Control	MCF-7, Jurkat.
Format	Liquid
Size	100 µl
Buffer	1×TBS (pH7.4), 1% BSA, 40% Glycerol.
Preservative	0.05% Sodium Azide
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw

BACKGROUND

Introduction

Acidic and basic fibroblast growth factors (FGFs) are members of a family of multifunctional polypeptide growth factors that stimulate proliferation of cells of mesenchymal, epithelial and neuroectodermal origin. Like other growth factors, FGFs act by binding and activating specific cell surface receptors. These include the Flg receptor or FGFR-1, the Bek receptor (or FGFR-2), FGFR-3, FGFR-4, FGFR-5 and FGFR-6. These receptors usually contain an extracellular ligand-binding region containing three immunoglobulin-like domains, a transmembrane domain and a cytoplasmic tyrosine kinase domain. The gene encoding human Bek (also designated K-sam) maps to chromosome 10q26.13 and is alternatively spliced to produce several isoforms. Heterogeneous mutations in Bek are associated with a range of craniosynostosis syndromes including Pfeiffer syndrome, Crouzon syndrome, Jackson-Weiss syndrome and Apert syndrome..

Keywords

bacteria-expressed kinase;BBDS;BEK;BEK fibroblast growth factor receptor;BFR1;CD332;CD332 antigen;CEK3;CFD1;Craniofacial dysostosis 1;ECT1;FGF receptor;FGFR 2;FGFR-2;Fgfr2;FGFR2_HUMAN;Fibroblast growth factor receptor 2;Hydroxyaryl protein kinase;Jackson Weiss syndrome;JWS;K SAM;K-sam;Keratinocyte growth factor receptor 2;Keratinocyte growth factor receptor;KGFR;KSAM;protein tyrosine kinase, receptor like 14;soluble FGFR4 variant 4;TK14;TK25 antibody
