



# Chimeric Rabbit Anti-Mouse CD147 monoclonal antibody, clone OX114 (CABT- L949R)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Chimeric Rabbit Anti-Mouse CD147 monoclonal antibody, clone OX114
<b>Specificity</b>	This clone binds to both the 2- and 3-domain isoforms of CD147.
<b>Immunogen</b>	Recombinant mCD147d0-2rCD4d3+4.
<b>Isotype</b>	IgG, κ
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Clone</b>	OX114
<b>Purification</b>	Protein A affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	FC, IHC, WB
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	200 µg, 1 mg
<b>Buffer</b>	PBS with 0.02% Proclin 300.
<b>Preservative</b>	0.02% Proclin 300.

<b>Storage</b>	Store at -20°C.
----------------	-----------------

<b>Ship</b>	Wet ice
-------------	---------

---

## BACKGROUND

<b>Keywords</b>	Basigin;basic immunoglobulin superfamily;HT7 antigen;membrane glycoprotein gp42;gp42;gp-42;OX114;MRC OX-114;MRC OX114
-----------------	---

---

## GENE INFORMATION

<b>Official Symbol</b>	basigin
------------------------	---------

<b>Synonyms</b>	BSG; basigin; HT-7; CD147; EMMPRIN; AI115436; AI325119; gp 42; HT7 antigen; neurothelin; 5A11/Basigin; membrane glycoprotein gp42; basic immunoglobulin superfamily;
-----------------	--