



Rabbit Anti-Human BTNL2 monoclonal antibody, clone FQS23363 (CABT-L962)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	BTNL2
Immunogen	Recombinant fragment within Human BTNL2 aa 200-400.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	FQS23363
Conjugate	Unconjugated
Applications	WB
Positive Control	Human spleen, tonsil, small intestine tissue lysates, K562 cell lysate.
Format	Liquid.
Size	100 µl
Buffer	59% PBS, 40% Glycerol, 0.05% BSA
Preservative	0.01% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Ship	Shipped at 4°C.

BACKGROUND

Introduction	Butyrophilin-like 2 (BTNL2; also BTL-II) is a predicted 60-65 kDa member of the BTN/MOG family, Ig-superfamily of proteins.
Keywords	BTNL2;butyrophilin-like 2 (MHC class II associated);SS2;BTN7;BTL-II;HSBLMHC1;butyrophilin-like protein 2

GENE INFORMATION

Entrez Gene ID	56244
UniProt ID	Q9UIR0