



## Mouse Anti-Human GDF-15 monoclonal antibody, clone 21B4 (CABT-L975M)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Specificity</b>	This antibody was selected for its ability to detect human GDF-15.
<b>Target</b>	Human GDF15
<b>Immunogen</b>	Human recombinant GDF-15
<b>Isotype</b>	IgG2
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	21B4
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, FC, Blocking, Neut, CyTOF-ready
<b>Reconstitution</b>	Lyophilized from a 0.2 um filtered solution in PBS.
<b>Format</b>	Lyophilized
<b>Size</b>	25 µg, 100 µg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C for up to 1 month. For longer storage, aliquot and store at -20°C.

Ship

Wet ice

---