



# Mouse Anti-Chicken CCL5 Monoclonal Antibody, clone 23I22 (CABT-NS1664)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	CCL5
Immunogen	Yeast expressed chicken CCL5
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Chicken
Clone	23I22
Purification	Protein G
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1µg/mL Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Size	100 µg
Buffer	PBS, 0.05% (w/v) Sodium Azide
Preservative	0.05% (w/v) Sodium Azide

**Storage**

-20°C

---

## **BACKGROUND**

**Keywords**

CCL5; chemokine (C-C motif) ligand 5; SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E;  
SIS-delta; C-C motif chemokine 5

---