



Chicken Anti-Human C-reactive Protein/CRP Polyclonal Antibody (CABT-NS1684)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CRP
Immunogen	Whole human C-Reactive Protein Mw= 23 KDa
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA, SPR Tests positively in indirect ELISA assay towards recombinant MRP 14 antigen when discriminated against a negative serum control. Specificity and stability validated by SPR. Apparent kD value estimated to be within a low nM (nanomolar) range. Uniformity between different production batches confirmed by SPR.
Format	Liquid, Purified
Size	1 ea
Buffer	PBS
Preservative	0.02% sodium azide
Storage	Store at 4°C short term. Longer term storage requires aliquoting and freezing at -20°C. Avoid freeze thaw cycles.

BACKGROUND

Keywords

CRP; C-reactive protein, pentraxin-related; PTX1; C-reactive protein; pentraxin 1;
