



Mouse Anti-BLV gp51 D-D' monoclonal antibody, clone CMW3 (CABT-RM001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	gp51 D-D'
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	BLV
Clone	CMW3
Conjugate	Unconjugated
Applications	RIPA, WB, FC, ELISA, IAC
Epitope	It recognizes the G epitope
Format	Liquid, Ascites
Concentration	1.0 mg/ml
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C. Do not freeze!

BACKGROUND

Introduction	51kDa surface glycoprotein
Keywords	Bovine leukemia virus; BLV; Bovine leukemia virus glycoprotein 51kDa; BLV glycoprotein 51kDa; Bovine leukemia virus gp51; BLV gp51