



Mouse Anti-BHV1/IBR gB-gI monoclonal antibody, clone E0F8 (CABT-RM002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	gB-gI
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	BHV1, BHV5
Clone	E0F8
Conjugate	Unconjugated
Applications	ELISA, Neut, IF
Format	Liquid, Ascites
Concentration	0.19 mg/ml
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	0.09% sodium azide
Storage	Store at 2-8°C. Do not freeze!

BACKGROUND

Introduction	Bovine herpesvirus 1 (BHV-1) is a virus of the family Herpesviridae and the subfamily Alphaherpesvirinae, known to cause several diseases worldwide in cattle, including
--------------	--

rhinotracheitis, vaginitis, balanoposthitis, abortion, conjunctivitis, and enteritis.

Keywords

Bovine Herpesvirus Type 1; BHV1; Bovine Herpesvirus Type 1 (IBR) glycoprotein B (GI); BHV1 glycoprotein B (GI); Bovine Herpesvirus Type 1 gB-gl; BHV1 gB-gl
