



Mouse Anti-CAEV 28 kDa capsid protein monoclonal antibody, clone DBFQ6B2 (CABT- RM004)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	28 kDa capsid protein
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	CAEV, MVV, OPPV
Clone	DBFQ6B2
Conjugate	Unconjugated
Applications	IP, IF
Format	Liquid, Ascites
Concentration	1.0 mg/ml
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	0.09% sodium azide
Storage	Store at 2-8°C. Do not freeze!

BACKGROUND

Introduction

CAEV 28 kDa capsid protein

Keywords

Caprine Arthritis Encephalitis Virus; CAEV; Caprine Arthritis Encephalitis Virus 28 kDa capsid protein; CAEV 28 kDa capsid protein