



Mouse Anti-N. caninum gp65 monoclonal antibody, clone 6C7-36 (CABT-RM005)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	gp65
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	N. caninum
Clone	6C7-36
Conjugate	Unconjugated
Applications	ELISA, IF, WB, Dot, PIHC
Format	Liquid, Ascites
Concentration	1.0 mg/ml
Size	100 µg
Buffer	phosphate-buffered saline
Preservative	0.09% sodium azide
Storage	Store at 2-8°C. Do not freeze!

BACKGROUND

Introduction	65 kDa glycosylated protein; structural phage protein
--------------	---

Keywords

Neospora caninum; N. caninum; Neospora caninum 65 kDa glycosylated protein; N. caninum 65 kDa glycosylated protein; Neospora caninum gp65; N. caninum gp65
