



Rabbit Anti-*C. albicans* monoclonal antibody (CABT-RM007)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	<i>C. albicans</i>
Isotype	IgG1
Source/Host	Rabbit
Species Reactivity	<i>C. albicans</i> (10 strains), <i>C. parapsilosis</i> , <i>C. glabrata</i> , <i>C. tropicales</i> , <i>C. auris</i>
Clone	7516
Purification	protein A chromatography
Conjugate	Unconjugated
Applications	IF, IHC, ELISA
Format	Liquid
Concentration	100 µg/ml
Size	100 µg
Buffer	0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide
Storage	Store at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C.

BACKGROUND

Introduction *Candida albicans* is an opportunistic pathogenic yeast that is a common member of the human

gut flora. It does not proliferate outside the human body. It is detected in the gastrointestinal tract and mouth in 40–60% of healthy adults.

Keywords

C. albicans; Candida albicans; Candida
