



Rabbit Anti-Aspergillus spp. monoclonal antibody (CABT-RM011)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Aspergillus
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Aspergillus
Clone	6256
Purification	protein A chromatography
Conjugate	Unconjugated
Applications	IF, ELISA
Format	Liquid
Concentration	100 µg/ml
Size	100 µg
Buffer	0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide
Storage	Store at 2-8°C until ready for use. Avoid repeated freeze-thawing by storing multiple aliquots at -20°C.

BACKGROUND

Introduction	Reactive with a secreted hyphal antigen from all Aspergillus species tested.
---------------------	--

Keywords

Aspergillus; fungus; mold
