



# Goat Anti-CCoV polyclonal antibody (CABT-RM019)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	CCV
Target	CCV
Source/Host	Goat
Species Reactivity	CCV
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, IP, WB
Format	Liquid
Concentration	1.0 mg/mL
Size	1 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.1% sodium azide
Storage	+4°C (2-8°C allowed)

## BACKGROUND

Introduction	Goats were immunized with Canine Coronavirus (CCV) nucleocapsid protein. After multiple immunizations in Freund's adjuvant, serum was collected and the anti-CCV immunoaffinity
--------------	---

purified.

---

**Keywords**

Canine Coronavirus; Coronaviru; Canine; CCV; Coronaviridae

---