



Mouse Anti-BCoV monoclonal antibody, clone 6B5 (CABT-RM020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	BCV
Target	BCV
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	BCV
Clone	6B5
Purification	Chromatography on protein G Sepharose
Conjugate	Unconjugated
Applications	ELISA, HIT
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage	+4°C (2-8°C allowed)

BACKGROUND

Introduction	Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with bovine corona virus.
--------------	-------------------------------------------------------------------------------------------------------------------------------------------------

Keywords

Bovine coronavirus; Bovine; coronavirus; BCV; Coronaviridae
