



Chicken Anti-MHV68 ORF59 monoclonal antibody (CABT-RM032)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	MHV68
Target	MHV68 ORF59
Immunogen	2 peptides (CGKKTRGGNKASDSGT, CKRPPPKKDREPTTKRPKL)
Isotype	IgY
Source/Host	Chicken
Species Reactivity	MHV68
Clone	BPSG60
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	10 mg
Buffer	Phosphate buffered saline, pH 7.2
Preservative	0.075% sodium azide
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles.

BACKGROUND

Introduction	DNA polymerase processivity factor
Keywords	Murine Gammaherpesvirus-68 ORF59; MHV68 ORF59; Murine; Gammaherpesvirus; ORF