



Mouse Anti-FOB Monoclonal antibody, clone C4006N (CABT-RM285)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Recognizes Human Hemoglobin
Target	FOB
Immunogen	Native Human Hemoglobin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	C4006N
Purification	≥ 90% pure by SDS-PAGE. Protein A Chromatography
Conjugate	unconjugated
Applications	LFIA, ELISA
Format	Liquid
Concentration	4.70 mg/mL
Size	1 mg
Buffer	Phosphate Buffered Saline, pH 7.4
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C. DO NOT FREEZE!

BACKGROUND

Introduction	Fecal occult blood (FOB) refers to blood in the feces that is not visibly apparent (unlike other types of blood in stool such as melena or hematochezia). A fecal occult blood test (FOBT) checks for hidden (occult) blood in the stool (feces).
Keywords	Fecal Occult Blood; FOB; Hemoglobin