



Mouse Anti-MPV NP Monoclonal antibody, clone C3973N (CABT-RM297)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Specificity | Recognizes Metapneumovirus Nucleoprotein |
| Target | MPV NP |
| Immunogen | Human MPV Strain NCL-174 |
| Isotype | IgG1 |
| Source/Host | Mouse |
| Species Reactivity | MPV |
| Clone | C3973N |
| Purification | > 95% pure (reduced SDS-PAGE). Protein A Chromatography |
| Conjugate | unconjugated |
| Applications | ELISA, IF, IHC, WB |
| Format | Liquid |
| Concentration | 1.0 mg/mL |
| Size | 1 mg |
| Buffer | Phosphate Buffered Saline, pH 7.4 |
| Preservative | 0.02% Sodium Azide |
| Storage | Store at 2-8°C for < 2 weeks. Long term storage at -20°C. Aliquot to avoid freeze/thaw cycles. |

BACKGROUND

| | |
|--------------|--|
| Introduction | Metapneumovirus is a genus of viruses in the family Pneumoviridae. |
| Keywords | Metapneumovirus Nucleoprotein; MPV NP; MPV; Metapneumovirus |